

Pillow talk



REUSABLE & DISPOSABLE **PILLOWS**
FOR STANDARD & SPECIALTY USE

good things in small packages

The Medline Space Saver Program*

Many Medline pillows undergo a unique packaging process that compresses the pillows enough for a 50 percent reduction in carton size, resulting in smaller cartons for delivery and space-saving storage for our customers. But rest assured – when removed from their packaging, the vacuum-packed pillows quickly plump back up to their full size for maximum patient comfort.

Disposable pillows



Fully compressed



Fully expanded
(within seconds after opening)

Reusable pillows



Fully compressed



Partially expanded
(within two hours after opening)



Fully expanded
(within four hours after opening)



Prior product carton, 12 each (back right)

NEW "compressed" product carton, 12 each (front left)



Look for this special "logo" to identify pillows packaged using Medline's **Space Saver** technology.

* Patent pending

Reusable Pillows




Nylex Ultra

Smooth, high-quality nylon ticking

- ❖ Vacuum packed two per poly bag to save much-needed storage space (*new packaging results in a 50 percent reduction in carton size*)
- ❖ Fully expands within four hours
- ❖ Fluid-proof woven nylon ticking
- ❖ Antimicrobial
- ❖ Inherently flame retardant
- ❖ Stain resistant
- ❖ Wipeable
- ❖ 18 oz polyester fill for 20" x 26",
16 oz polyester fill for 18" x 24"



<i>Item No.</i>	<i>Description</i>	<i>Qty</i>
MDT219715	Tan, 20" x 26"	12/cs
MDT219716	Tan, 18" x 24"	20/cs



Ovation®

Great for the budget-conscious

- ❖ Vacuum packed two per poly bag to save much-needed storage space (*new packaging results in a 50 percent reduction in carton size*)
- ❖ Fully expands within four hours
- ❖ Soft, fluid-proof polypropylene ticking
- ❖ Hypo-allergenic ticking
- ❖ Antimicrobial
- ❖ Inherently flame retardant
- ❖ Stain resistant
- ❖ Wipeable
- ❖ 18 oz polyester fill for 20" x 26",
16 oz polyester fill for 18" x 24"



<i>Item No.</i>	<i>Description</i>	<i>Qty</i>
MDT219860	White, 20" x 26"	12/cs
MDT219861	White, 18" x 24"	20/cs
MDT219885	Blue, 20" x 26"	12/cs
MDT219886	Blue, 18" x 24"	20/cs


Reusable Pillows

Medsoft®

A soft and comfortable choice

- ❖ Mold-resistant vinyl ticking
 - ❖ Antimicrobial
 - ❖ Inherently flame retardant
 - ❖ Fluid-proof
 - ❖ Antistatic
 - ❖ Stain resistant
 - ❖ Wipeable
- ❖ 18 oz polyester fill for 20" x 26",
16 oz polyester fill for 18" x 24"

Item No.	Description	Qty
MDT219683	White, 20" x 26"	12/cs
MDT219684	White, 18" x 24"	20/cs



Nylex II

Same great nylon ticking in different options

- ❖ Fluid-proof woven nylon ticking
 - ❖ Antimicrobial
 - ❖ Inherently flame retardant
 - ❖ Stain resistant
 - ❖ Wipeable
- ❖ 18 oz polyester fill for 20" x 26",
16 oz polyester fill for 18" x 24"

Item No.	Description	Qty
MDT219715D	Tan, 20" x 26"	12/cs
MDT219715GH	Tan 20" x 26" w/grommet	12/cs
MDT219716D	Tan, 18" x 24"	20/cs
MDT219800D	Blue, 20" x 26"	12/cs
MDT219801D	Blue, 18" x 24"	20/cs
MDT219717D	Tan, 13" x 17" positioning	24/cs
MDT219718D	Tan, 9" x 12" positioning	12/cs



Metal grommet allows for convenient pillow hanging and storage



Reusable Pillows



Royal Diamond®

A rich, fluffy at-home feel

- ❖ Part of our exclusive *Feels Like Home™* line
- ❖ Soft 100 percent cotton ticking (330 thread count)
- ❖ Hypo-allergenic
- ❖ Machine washable
- ❖ 18 oz polyester fill for 20" x 26",
22 oz polyester fill for 20" x 30"

<i>Item No.</i>	<i>Description</i>	<i>Qty</i>
MDT219601	White, 20" x 26"	12/cs
MDT219600	White, 20" x 30"	10/cs




Stay-Fluff

Foam core for extra durability

- ❖ High-density foam core wrapped in 14 oz of 100 percent shape-maintaining garneted polyester
- ❖ Antimicrobial
- ❖ Fluid-proof woven nylon ticking
- ❖ Inherently flame retardant
- ❖ Stain resistant
- ❖ Wipeable

<i>Item No.</i>	<i>Description</i>	<i>Qty</i>
MDT219721D	Tan, 20" x 26"	12/cs



Foam core helps pillow maintain its shape


Reusable Pillows

Velvet Soft

A soft, launderable option

- ❖ Machine washable for easy maintenance
 - ❖ Antimicrobial
 - ❖ Fluid resistant
- ❖ 18 oz polyester fill for 20" x 26",
16 oz polyester fill for 18" x 24"

<i>Item No.</i>	<i>Description</i>	<i>Qty</i>
MDT219708D	White, 20" x 26"	12/cs
MDT219707	White, 18" x 24"	20/cs



Hyperbaric

Great for your hyperbaric requirements

- ❖ Suitable for any hyperbaric chamber
 - ❖ 100 percent cotton ticking in blue and white stripes
- ❖ Label-free with 100 percent cotton components for hyperbaric safety
- ❖ 30 oz 100 percent garnetted cotton fill

<i>Item No.</i>	<i>Description</i>	<i>Qty</i>
MDTALLCOTT2	Blue-and-white striped, 20" x 26"	12/cs



Disposable Pillows



Classic

Helps reduce the risk of cross-contamination

- ❖ Three individually bagged pillows per vacuum-packed master poly bag
- ❖ Fully expands within seconds
- ❖ Soft 100 percent olefin ticking
- ❖ Lightweight and durable
- ❖ Prohibits the growth of mildew
- ❖ Inherently flame-retardant



<i>Item No.</i>	<i>Description</i>	<i>Qty</i>
NON24392	White, 12" x 16", lt. wt.	24/cs
NON24391	White, 16" x 22", lt. wt.	12/cs
NON2439114*	White, 16" x 22", med. wt.	12/cs
NON24390	White, 18" x 24", med. wt.	12/cs
NON24397*	White, 19" x 26", med. wt.	12/cs
NON24393	White, 21" x 27", med. wt.	12/cs
NON2439322	White, 21" x 27", hvy. wt.	12/cs

* Not vacuum packed



ComfortMed®

Ideal for the ER and OR

- ❖ Our most popular option
- ❖ Three individually bagged pillows per vacuum-packed master poly bag
- ❖ Fully expands within seconds
- ❖ Soft 100 percent olefin white ticking with blue hem thread
- ❖ Lightweight and durable
- ❖ Prohibits the growth of mildew
- ❖ Inherently flame-retardant



<i>Item No.</i>	<i>Description</i>	<i>Qty</i>
PM1622-10	White, 16" x 22", lt. wt.	12/cs
PM1824-15	White, 18" x 24", med. wt.	12/cs
PM2127-18	White, 21" x 27", med. wt.	12/cs
PM2127-22	White, 21" x 27", hvy. wt.	12/cs

ordering information

STANDARD REUSABLE

	Item No.	Description	Qty
<i>Nylex Ultra</i>	MDT219715 SS	Tan, 20" x 26" (51cm x 66cm)	12/cs
	MDT219716 SS	Tan, 18" x 24" (46cm x 61cm)	20/cs
<i>Ovation®</i>	MDT219860 SS	White, 20" x 26" (51cm x 66cm)	12/cs
	MDT219861 SS	White, 18" x 24" (46cm x 61cm)	20/cs
	MDT219885 SS	Blue, 20" x 26" (51cm x 66cm)	12/cs
	MDT219886 SS	Blue, 18" x 24" (46cm x 61cm)	20/cs
<i>Medsoft®</i>	MDT219683	White, 20" x 26" (51cm x 66cm)	12/cs
	MDT219684	White, 18" x 24" (46cm x 61cm)	20/cs

SPECIALTY REUSABLE

<i>Nylex II</i>	MDT219715D	Tan, 20" x 26" (51cm x 66cm)	12/cs
	MDT219715GH	Tan, 20" x 26" (51cm x 66cm) w/grommet	12/cs
	MDT219716D	Tan, 18" x 24" (46cm x 61cm)	20/cs
	MDT219800D	Blue, 20" x 26" (51cm x 66cm)	12/cs
	MDT219801D	Blue, 18" x 24" (46cm x 61cm)	20/cs
	MDT219717D	Tan, 13" x 17" (33cm x 43cm) positioning	24/cs
	MDT219718D	Tan, 9" x 12" (23cm x 30cm) positioning	12/cs
<i>Royal Diamond®</i>	MDT219601	White, 20" x 26" (51cm x 66cm)	12/cs
	MDT219600	White, 20" x 30" (51cm x 76cm)	10/cs
<i>Stay-Fluff</i>	MDT219721D	Tan, 20" x 26" (51cm x 66cm)	12/cs
<i>Velvet Soft</i>	MDT219708D	White, 20" x 26" (51cm x 66cm)	12/cs
	MDT219707	White, 18" x 24" (46cm x 61cm)	20/cs
<i>Hyperbaric</i>	MDTALLCOTT2	Blue and White Stripe, 20" x 26" (51cm x 66cm)	12/cs

DISPOSIBLE

<i>Classic</i>	NON24392 SS	White, 12" x 16" (30cm x 41cm), lt. wt.	24/cs
	NON24391 SS	White, 16" x 22" (41cm x 56cm), lt. wt.	12/cs
	NON2439114	White, 16" x 22" (41cm x 56cm), med. wt.	12/cs
	NON24390 SS	White, 18" x 24" (46cm x 61cm), med. wt.	12/cs
	NON24397	White, 19" x 26" (48cm x 66cm), med. wt.	12/cs
	NON24393 SS	White, 21" x 27" (53cm x 69cm), med. wt.	12/cs
	NON2439322 SS	White, 21" x 27" (53cm x 69cm), hvy. wt.	12/cs
<i>ComfortMed®</i>	PM1622-10 SS	White, 16" x 22" (41cm x 56cm), lt. wt.	12/cs
	PM1824-15 SS	White, 18" x 24" (46cm x 61cm), med. wt.	12/cs
	PM2127-18 SS	White, 21" x 27" (53cm x 69cm), med. wt.	12/cs
	PM2127-22 SS	White, 21" x 27" (53cm x 69cm), hvy. wt.	12/cs

SS Space Saver packed

1-800-MEDLINE (633-5463)
www.medline.com

Medline Canada 1-800-396-6996
www.medline.ca | canada@medline.com

Medline México 01-800-831-0898
www.medlinemexico.com | mexico@medline.com

©2009 Medline Industries, Inc.
ComfortMed, Medline, Medsoft and Ovation are registered trademarks of Medline Industries, Inc. Royal Diamond is a registered trademark of Carpenter Company.
We reserve the right to correct any errors that may occur within this brochure.

Some products may not be available for sale in Canada or Mexico. Please contact your Account Representative for additional information.

Medline Industries, Inc.



One Medline Place | Mundelein, IL 60060

MKT209049 / LT423 / 3FE / 10M